

Anti-CD154 antibody [MK13A4] {FITC} (STJ16100541)

STJ16100541

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD154 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|------------------------------|
| Clonality | Monoclonal |
| Clone ID | MK13A4 |
| Concentration | |
| Conjugation | FITC |
| Purification | |
| Dilution Range | |
| Formulation | TBS and 1% BSA |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |