

Anti-CD236 antibody [BRIC100] (STJ16100543)

STJ16100543

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD236 is suitable for use in Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	WB, IHC-F, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BRIC100
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TRIS buffered glycine
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence