

Anti-CD173 antibody [BRIC231] (STJ16100544)

STJ16100544

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD173 is suitable for use in Immunofluorescence and Flow Cytometry research applications.
Applications	IF, FC
Host/Source	Mouse
Reactivity	Human, Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BRIC231
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TRIS buffered glycine
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	