

Anti-CD88 antibody [P12/1] (STJ16100547)

STJ16100547

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD88 is suitable for use in Western Blot, Immunohistochemistry, Flow Cytometry and Immunohistochemistry research applications.
Applications	WB, IHC-P, FC, IHC-F
Host/Source	Mouse
Reactivity	Human, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	P12/1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence