

Anti-CD229 antibody [HLy9.1.25] {PE} (STJ16100551)

STJ16100551

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD229 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HLy9.1.25
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	