

Anti-CD101 antibody [BB27] (STJ16100560)

STJ16100560

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD101 is suitable for use in Immunoprecipitation, Flow Cytometry and Immunohistochemistry research applications.
Applications	IP, FC, IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BB27
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	9398
Gene Symbol	CD101
Uniprot ID	IGSF2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	