

Anti-CD101 antibody [BB27] {FITC} (STJ16100561)

STJ16100561

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD101 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BB27
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	9398
Gene Symbol	CD101
Uniprot ID	IGSF2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	