

Anti-CD226 antibody [DX11] {Biotin} (STJ16100569)

STJ16100569

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD226 is suitable for use in Immunoprecipitation and Flow Cytometry research applications.
Applications	IP, FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DX11
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	10666
Gene Symbol	CD226
Uniprot ID	CD226_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	