

Anti-CD226 antibody [DX11] {PE} (STJ16100571)

STJ16100571

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD226 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | DX11 |
| Concentration | |
| Conjugation | PE |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.09% sodium azide and 1% BSA, 5% sucrose |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 10666 |
| Gene Symbol | CD226 |
| Uniprot ID | CD226_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |