

## Mouse Anti-Human IgA antibody {FITC} (STJ16100581)

STJ16100581

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	FITC-conjugated mouse monoclonal anti-Human IgA secondary antibody. For use in most research applications.
<b>Applications</b>	IHC, ELISA, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	BFG-267
<b>Concentration</b>	1.2 mg/ml
<b>Conjugation</b>	FITC
<b>Purification</b>	
<b>Dilution Range</b>	IHC ELISA 1:1000 WB
<b>Formulation</b>	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	973
<b>Gene Symbol</b>	CD79A
<b>Uniprot ID</b>	CD79A_HUMAN
<b>Immunogen</b>	IgA
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	