

Anti-Tumor Necrosis Factor Alpha antibody [52B83] (STJ16100583)

STJ16100583

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tumor Necrosis factor Alpha is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA and Western Blot research applications.
Applications	IHC-F, FC, ELISA, WB
Host/Source	Mouse
Reactivity	Human, Guinea Pig, Mouse, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	52B83
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence