

Anti-Immunoglobulin IgA antibody [MH14/1] (STJ16100584)

STJ16100584

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Immunoglobulin IgA is suitable for use in Immunohistochemistry, ELISA and Immunofluorescence research applications. |
| Applications | IHC-F, IHC-P, ELISA, IF |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Monoclonal |
| Clone ID | MH14/1 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | IHC-F IHC-P (1:5-1:10) ELISA IF |
| Formulation | Contains 0.09% sodium azide |
| Isotype | IgG1 |
| Storage | Store at 2-8°C upon receipt. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |