

Mouse Anti-Human IgG antibody {HRP} (STJ16100589)

STJ16100589

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated mouse monoclonal anti-Human IgG secondary antibody. For use in most research applications.
Applications	WB, ELISA, IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EFE-565
Concentration	1.2 mg/ml
Conjugation	HRP
Purification	
Dilution Range	WB ELISA IHC
Formulation	0.05 M Sodium Phosphate, 0.1 M NaCl, 50% Glycerol, 5 mg/ml BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	IgG
Immunogen Region	
Specificity	
Immunogen Sequence	