

Anti-Interleukin 6 antibody [MP5-32C11] (STJ16100594)

STJ16100594

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin 6 is suitable for use in ELISA and Immunoprecipitation. Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MP5-32C11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
Uniprot ID	IL6_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	