

Anti-Urokinase antibody [PGM2001] (STJ16100595)

STJ16100595

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Urokinase is suitable for use in Western Blot, ELISA and Immunohistochemistry research applications.
Applications	WB, ELISA, IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PGM2001
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB (1:100) ELISA (1:5000) IHC-F (1:30)
Formulation	Purified and lyophilized
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence