

Anti-LDL Apolipoprotein-E antibody [3D12] (STJ16100606)

STJ16100606

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-LDL Apolipoprotein-E is suitable for use in Immunohistochemistry, ELISA and Western Blot research applications.
Applications	IHC, ELISA, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3D12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB ELISA (x1000) IHC paraffin sections (x100)
Formulation	Purified and lyophilized
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence