

Mouse Anti-Tag Biotin antibody {Biotin} (STJ16100612)

STJ16100612

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Secondary antibodies |
| Short Description | Biotin-conjugated mouse monoclonal anti-Tag Biotin secondary antibody. For use in most research applications. |
| Applications | IHC, ELISA, FC |
| Host/Source | Mouse |
| Reactivity | Tag |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | Bio-8 |
| Concentration | 100 ug/ml |
| Conjugation | Biotin |
| Purification | |
| Dilution Range | IHC ELISA FC |
| Formulation | Tissue culture supernatant with Sodium azide |
| Isotype | IgG3 |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |