

Mouse Anti-Tag Streptavidin antibody {Unconjugated} (STJ16100613)

STJ16100613

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Unconjugated mouse monoclonal anti-Tag Streptavidin secondary antibody. For use in most research applications.
Applications	IHC, ELISA
Host/Source	Mouse
Reactivity	Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Strep-10
Concentration	200 ug/ml
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:20 ELISA
Formulation	Supernatant with 1% BSA and 0.1% Sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence