

## Mouse Anti-Chicken IgY antibody {Biotin} (STJ16100615)

STJ16100615

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Biotin-conjugated mouse monoclonal anti-Chicken IgY secondary antibody. For use in most research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Chicken

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	7C2
<b>Concentration</b>	50 ug/ml
<b>Conjugation</b>	Biotin
<b>Purification</b>	
<b>Dilution Range</b>	ELISA, WB
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% Sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen Sequence</b>