

Rat Anti-Mouse IgM antibody {HRP} (STJ16100645)

STJ16100645

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated rat monoclonal anti-Mouse IgM secondary antibody. For use in most research applications.
Applications	IHC-F, ELISA
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LO-MM-9
Concentration	1 mg/ml
Conjugation	HRP
Purification	
Dilution Range	IHC-fr ELISA
Formulation	PBS with 0.1% Sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	NA
Uniprot ID	IGM_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	