

Anti-Interleukin-1 receptor/P80 clone Reg21 antibody [Reg21] (STJ16100655)

STJ16100655

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Interleukin-1 receptor/P80 clone Reg21 is suitable for use in ELISA research applications. |
| Applications | ELISA |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | Reg21 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.1% BSA and 0.02% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |