

Anti-IFN-gamma antibody [F14] (STJ16100661)

STJ16100661

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IFN-gamma is suitable for use in Immunohistochemistry, ELISA and Immunoprecipitation research applications.
Applications	IHC-F, ELISA, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	F14
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3458
Gene Symbol	IFNG
Uniprot ID	IFNG_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	