

## Anti-IFN-Alpha antibody [N27] (STJ16100665)

STJ16100665

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-IFN-Alpha is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	ELISA, IHC-F, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	N27
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	ELISA WB (1-2 Mu g/ml) IHC frozen sections (1-2 Mu g/mlfor 30 min at RT).
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

