

## Anti-CD90 antibody [F15-42-1] (STJ16100673)

STJ16100673

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD90 is suitable for use in Immunofluorescence, Immunoprecipitation, Flow Cytometry and Immunohistochemistry research applications.
<b>Applications</b>	IF, IP, FC, IHC-F
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Cynomolgus, Monkey

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	F15-42-1
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>