

Anti-CD90 antibody [F15-42-1] {FITC} (STJ16100674)

STJ16100674

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD90 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human, Cynomolgus, Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	F15-42-1
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	