

Anti-CD90 antibody [F15-42-1] {PE} (STJ16100675)
STJ16100675**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD90 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human, Cynomolgus, Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	F15-42-1
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA, 5% sucrose
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	