

Anti-Endothelial Cell Antigen-1 antibody [HIS-52] (STJ16100679)

STJ16100679

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Endothelial Cell antigen-1 is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-F, FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HIS-52
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Supernatant with FCS, with azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence