

## Anti-E-Selectin CD62-E antibody [ENA1] (STJ16100682)

STJ16100682

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-E-Selectin CD62-E is suitable for use in Immunohistochemistry, ELISA and Immunoprecipitation research applications.
<b>Applications</b>	IHC-F, ELISA, IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Baboon, Chimpanzee, Cynomolgus, Monkey, Rhesus Monkey

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	ENA1
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC 1:10 WB 1:10 ELISA 1:10 IP 1:10
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	