

Anti-E-Selectin CD62-E antibody [ENA2] (STJ16100683)

STJ16100683

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-E-Selectin CD62-E is suitable for use in Immunohistochemistry, ELISA and Immunofluorescence research applications.
Applications	IHC-F, ELISA, IF
Host/Source	Mouse
Reactivity	Human, Baboon, Chimpanzee, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ENA2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence