

Anti-VCAM-1 antibody [1G11B1] {Biotin} (STJ16100685)

STJ16100685

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-VCAM-1 is suitable for use in research applications.
Applications	
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1G11B1
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7412
Gene Symbol	VCAM1
Uniprot ID	VCAM1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	