

Anti-NHERF antibody [EBP-10] (STJ16100695)

STJ16100695

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-NHERF is suitable for use in Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	IP, WB, IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBP-10
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1k
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	9368
Gene Symbol	SLC9A3R1
Uniprot ID	NHRF1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	