

Anti-MUC5AC antibody [45M1] (STJ16100709)

STJ16100709

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MUC5AC is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications.
Applications	ELISA, IHC-F, IHC-P, WB
Host/Source	Mouse
Reactivity	Human, Mouse, Rat, Rabbit, Monkey, Chicken, Cow, Cat, Pig, Hedgehog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	45M1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4586
Gene Symbol	MUC5AC
Uniprot ID	MUC5A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

