

Anti-TIMP-1 antibody [2A5] (STJ16100714)

STJ16100714

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TIMP-1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2A5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7076
Gene Symbol	TIMP1
Uniprot ID	TIMP1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	