

Anti-TIMP-2 antibody [3A4] (STJ16100715)

STJ16100715

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TIMP-2 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P, IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3A4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7077
Gene Symbol	TIMP2
Uniprot ID	TIMP2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	