

## Anti-CD45 antibody [ANK74] (STJ16100717)

STJ16100717

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD45 is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
<b>Applications</b>	IHC-F, FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	ANK74
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Culture supernatant with 0.05% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	5788
<b>Gene Symbol</b>	<a href="#">PTPRC</a>
<b>Uniprot ID</b>	<a href="#">PTPRC_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	