

Anti-T-cell Receptor antibody [R73-2b] (STJ16100718)

STJ16100718

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-T-cell Receptor is suitable for use in Immunohistochemistry and Flow Cytometry research applications. |
| Applications | IHC-F, FC |
| Host/Source | Mouse |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | R73-2b |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | Culture supernatant with 0.05% sodium azide |
| Isotype | IgG2b |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |