

Anti-T-cell Receptor & Epithelial Tumours antibody [R73xCC2] (STJ16100721)

STJ16100721

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell Receptor & Epithelial Tumours is suitable for use in research applications.
Applications	
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R73xCC2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	