

Anti-Influenza A Nucleoprotein antibody [AA5H] (STJ16100728)

STJ16100728

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Influenza A Nucleoprotein is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB, IHC-P, IF
Host/Source	Mouse
Reactivity	Viral

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AA5H
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence