

Anti-CD6 antibody [MEM-98] (STJ16100729)

STJ16100729

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD6 is suitable for use in Flow Cytometry, Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	FC, IP, WB, IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-98
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	923
Gene Symbol	CD6
Uniprot ID	CD6_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	