

Anti-Melanoma Associated Antigen antibody [NKI/M6] (STJ16100733)

STJ16100733

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Melanoma Associated antigen is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NKI/M6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F (1:10)
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	