

Anti-Vascular Endothelial Cadherin antibody [BV9] (STJ16100741)

STJ16100741

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Vascular Endothelial Cadherin is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F, FC, ELISA, IF, IP, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BV9
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1003
Gene Symbol	CDH5
Uniprot ID	CADH5_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	