

Anti-CD21 antibody [4C5] (STJ16100747)

STJ16100747

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD21 is suitable for use in Immunohistochemistry, Immunofluorescence and Flow Cytometry research applications.
Applications	IHC-F, IF, FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4C5
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F (1:20-1:50) IF FC (1:50)
Formulation	serum with 0.1% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	