

Anti-Granzyme B antibody [GrB-7] (STJ16100752)

STJ16100752

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Granzyme B is suitable for use in Immunocytochemistry, ELISA and Immunohistochemistry research applications. |
| Applications | ICC, ELISA, IHC-P |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Monoclonal |
| Clone ID | GrB-7 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | ICC 1:20 |
| Formulation | Serum free tissue culture supernatant with 0.1% BSA and 0.09% sodium azide |
| Isotype | IgG2a |
| Storage | Store at 2-8°C upon receipt. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|----------------------------|
| Gene ID | 3002 |
| Gene Symbol | GZMB |
| Uniprot ID | GRAB_HUMAN |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |