

## Anti-MCP-1 antibody [MNA.1] (STJ16100753)

STJ16100753

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-MCP-1 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC-F, FC, ELISA, IP, IHC-P, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Pig

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MNA.1
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">6347</a>
<b>Gene Symbol</b>	<a href="#">CCL2</a>
<b>Uniprot ID</b>	<a href="#">CCL2_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	