

Anti-EBV-EA P55 antibody [1108-1] (STJ16100758)

STJ16100758

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-EBV-EA P55 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, IHC-F, IF, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1108-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

