

Anti-p53 antibody [Pab122] (STJ16100763)

STJ16100763

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-p53 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IF
Host/Source	Mouse
Reactivity	Human, Mouse, Rat, Dog, Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Pab122
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7157
Gene Symbol	TP53
Uniprot ID	P53_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	

