

Anti-LBP antibody [M330-19] (STJ16100765)

STJ16100765

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-LBP is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M330-19
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3929
Gene Symbol	LBP
Uniprot ID	LBP_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	