

Anti-Rabies Virus antibody [RAB-50] (STJ16100770)

STJ16100770

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Rabies Virus is suitable for use in ELISA, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	ELISA, IHC-F, IF, WB
Host/Source	Mouse
Reactivity	Rabies Virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RAB-50
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

