

Anti-PR3 antibody [PRG-2] (STJ16100777)

STJ16100777

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-PR3 is suitable for use in Flow Cytometry, ELISA and Western Blot research applications.
Applications	FC, ELISA, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PRG-2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	5657
Gene Symbol	PRTN3
Uniprot ID	PRTN3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	