

Anti-RBP antibody [G4E4] (STJ16100778)

STJ16100778

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-RBP is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications. |
| Applications | ELISA, IHC-F, IHC-P, WB |
| Host/Source | Mouse |
| Reactivity | Chimpanzee, Goat, Human, Monkey, Mouse, Rabbit, Rat, Rhesus |

PRODUCT PROPERTIES

| | |
|----------------------------|------------------------------|
| Clonality | Monoclonal |
| Clone ID | G4E4 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.02% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |

